

# Roller Mixer/Neutral Standard Tube Roller MX-T6-S+

## Application

innoLAB laboratory tube roller provide a gentle and highly efficient rolling motions, deal for mixing blood samples, viscous substances and liquid-solid suspensions. Used in a variety of applications, such as the prevention blood coagulation, immune precipitation, etc.



| Specification                                | MX-T6-S+            |
|--|---------------------|
| Amplitude                                    | 24mm                |
| Motor type                                   | Brushless DC motor  |
| Rocking motion                               | Rocking and rolling |
| Max. load capacity                           | 4Kg                 |
| Number of roller                             | 6                   |
| Roller size [length]                         | 280mm               |
| Speed range                                  | 0 - 70rpm           |
| Speed display                                | Scale               |
| Operation mode                               | Continuous          |
| Voltage, Frequency                           | 100 - 240V, 50/60Hz |
| Power  | 25W                 |
| Dimension [D x W x H]                        | 485 x 250 x 100mm   |
| Weight                                       | 4.5Kg               |
| Permissible ambient temperature and humidity | 5 - 40°C, 80%RH     |
| Protection class                             | IP21                |

## Features

- Adjustable speed range from 0 to 70rpm
- Six roller design for gentle and efficient rolling, and is easy to clean

Manufactured by:

**DLAB SCIENTIFIC CO.,LTD**

31 Yu An Road, Beijing Airport Economic Core Zone,  
Shunyi District, Beijing 101318, China

Distributed by:

**PT INNOVASI DIAGNOSTIKA**

Kawasan Industri Jababeka III, Jl. Tekno 1 Blok C2 Unit C,  
Cikarang Utara 17530, Indonesia